

L1A FITZMILL® / QUADRO® COMIL®

IDEAL FOR LABORATORY & SMALL CAPACITY APPLICATIONS

- Incorporates the same process advantages of larger Fitzmill® Comminutors and Quadro® Conical Screen Mills.
- Easily converts from Fitzmill® Communitor head to Quadro® Comil® U5 head, with minimal effort.
- Easily integrates into glove box isolator and other containment systems.
- Can process samples as small as 5 grams with minimal loss of fines.
- Can process wet and/or dry products.
- Underdriven design provides maximum inlet opening for efficient feeding.
- Results are directly scalable to large FitzMill® Communitors and Quadro® Conical Screen Mills.
- The variable speed rotor/impeller provides the equivalent tip speeds of larger units.
- A wide variety of screens and rotors/impellers are available for process optimization.

ADDITIONAL OPTIONS INCLUDE:

- Base to support machine.
- Various size containers.
- Transformer for 110V operation.



832 Industrial Drive
Elmhurst, IL 60126 USA
Phone: 630-530-3333
Fax: 630-530-0832
Email: fitzinfo@idexcorp.com
Web: www.fitzmill.com



QUADRO
Leading Process Equipment Innovation

613 Colby Drive
Waterloo, ON, CA N2V 1A1
Phone: 519-884-9660
Fax: 519-884-0253
Email: quadrosales@idexcorp.com
Web: www.quadro.com



1925 West Field Court, Suite 200, Lake Forest, IL 60045-4824

L1A-QCU5/12